3D Direct Laser Writing



3D Direct Laser Writing is a tool to fabricate 3D freeform structures down sub μ m. It is based on *Two-Photon Lithography* but beyond that 2D and 2.5D structures with nano dimensions are also possible. This system uses a nonlinear two-photon absorption process to modify, *e.g.* polymerize, a photosensitive medium at a specific point in the resist. By scanning the photoresist with a stage over this point a 3D-structure with dimensions in the submicron scale or greater can be written.

Contact

See KNMF website or contact the KNMF User Office.

Equipment

- Nanoscribe Photonic Professional GT
 - Several sample holders
 - Galvoscan unit
 - Hybrid stage for large accurate travel distances
- Critical point dryer (supercritical CO₂)
- Coming soon: UV flood exposure (shell writing mode for large structures)

Features

- Resolution: 3D:200 nm lateral, 750 nm normal 2D: 180 nm
- Feature Structure size: max. 600 x 600 x 3700 μm³ depending on filling factor
- Writing modes: piezo (100x DIP), galvoscan (63x GT, 25x GT)
- Corresponding writing fields: 300 x 300 x 300 μm³; 140 x 140 μm²; 280 x 280 μm²
- Writing times: piezo slow, galvoscan fast
- Larger areas have to be stitched
- Accessible writing area: 100 x 100 mm² where structures could be placed


Limitations/constraints


- The realizable structure size depends on the structural stability of the design
- The best results can be reached by using the IPresists

Materials

• Negative / positive / experimental resists Details and substrates see next page

Sample structures





Karlsruhe Nano Micro Facility (KNMF) Karlsruhe Institute of Technology (KIT) Hermann-von-Helmholtz-Platz 1 76344 Eggenstein-Leopoldshafen Germany

User Office

 Phone:
 +49 (721) 608-23123

 Fax:
 +49 (721) 608-26273


 Email:
 knmf-useroffice@kit.edu

 Web:
 www.knmf.kit.edu

3D Direct Laser Writing



Writing principle



Sketch of the writing of a line in 3D-space inside a resist layer. The inset shows the modification in the voxel (blue) at the focal position. Only in the voxel two-photon absorption occurs.

Writing modes



- structure height
- Fixed thin substrate
- any resist usable



- Higher AR of voxel
- Bigger voxelBigger processing
- effort
 - Only solid resist
 - no transparent substrate neededHigh working


distance



- Not every resist usable
- Special glass substratesSmall working distance
- structure height
- Very small voxel



- Special glass substrates
- no vertical lines possible
- 100x faster than other writing modes





Resists

- Negative Resists: IP-L, IP-G, IP-Dip, IP-S, and similar resists that are photosensitive at a wavelength of 380nm
- Positive Resist: AZ9260 (in preparation)
- Experimental Resists are possible, but only the Air mode objective or oil-immersion Objective are applied. Dip-in techniques can only be used with proven compatibility.

Karlsruhe Nano Micro Facility (KNMF) Karlsruhe Institute of Technology (KIT) Hermann-von-Helmholtz-Platz 1 76344 Eggenstein-Leopoldshafen Germany

User Office

 Phone:
 +49 (721) 608-23123

 Fax:
 +49 (721) 608-26273

 Email:
 knmf-useroffice@kit.edu

 Web:
 www.knmf.kit.edu

3D Direct Laser Writing



Substrates

- 25 x 25 x 0.7 mm Glass, glass covered with ITO, cover slides 22 x 22 x 0.17 mm
- Si-wafer 4" (100 mm)
- Si/SiO2-wafer 4" (100 mm)
- Metallized Si-Wafer (Cr/Au)
- Other substrates have to be provided by the user.

Data

- Stl-format
- CAD data can be exported in stl-format:
- Sometimes CAD programs export erroneous stl-files. Solid works and similar programs used in mechanical engineering produce correct stl.files.
- To avoid errors within stl-files the following rules can help:
 - Use correct units during contruction (μ m), otherwise resolution could be bad
 - Use volumes, no areas, these objects have to be closed
 - Avoid Boolean operations
 - Avoid duplexes (use snaps)
 - Think about plane orientation. Normals are used to define interior or exterior. Therefore contruct plane consistently.
 - Avoid additional structures not necessary (not written) to your part. Correction of the stl-file is nearly impossible
 - In case of periodic structures a single unit cell is sufficient.

Karlsruhe Nano Micro Facility (KNMF) Karlsruhe Institute of Technology (KIT) Hermann-von-Helmholtz-Platz 1 76344 Eggenstein-Leopoldshafen Germany

User Office

 Phone:
 +49 (721) 608-23123

 Fax:
 +49 (721) 608-26273

 Email:
 knmf-useroffice@kit.edu

 Web:
 www.knmf.kit.edu